

1 KSK RUBBER DAM CLAMP

Autoclave 132°C

We have 3 surface types of rubber dam clamp.


Normal : Most common type.

Matt : Lusterless type.

Black : Lusterless type.

These new clamp "Black" and "MATT" protected from eye fatigue of deintists by reflecting strong operation light.

The surface color of "Black" is not paint but oxide layer color, so it is safe for human body.



2 TIGER CLAMP

Autoclave 132°C

"TIGER CLAMP" has teeth on their beak, and strong holding power like a tiger.

Bicuspid : #1T, #1AT, #2T, #2AT

Molar : #14T, #56T

Anterior : #9T

#1T : For roots deep festooned beaks.

#1AT : For cervical cavities


#2T : For lower bicuspids with flat beaks.

#2AT : Similar to #2T but larger beaks.

#14T : For half erupted molars.

#56T : For lower molars with festooned beaks.

#9T : Most common anterior clamp.



3 DECIDUOUS CLAMP

Autoclave 132°C

Compact and anatomical shape enables surely clamping.

#54 : For upper right D

#55 : For upper right E

#64 : For upper left D

#65 : For upper left E

#74 : For lower left D

#75 : For lower left E

#85 : For lower left E

#84 : For lower right D

#54



#55



#64



#65



#74



#75



#84




#85



4 DECIDUOUS CLAMP BOARD

Autoclave 132°C

All in one set of deciduous clamps.



5 MOLAR CLAMPS

Autoclave 132°C

We have both "Ivory type" and "White type" clamps.

#26: For upper molars especially for erupted tooth.

#28: For upper or lower molars. Very wide fringe.

#30: For beaked molars.

#31: For beaked molars.

#51: For beaked molars with both side beaks.

#138: For third molars.

#139: For third molars.

#200: For lower molars with wings and broad fringe.

#201: Similar to #200 with festooned beaks.

#202: For large molar, similar to #200, but with larger beaks.

#203: For small lower right and left molars.

#204: For small lower right and left molars.

#205: For upper molar with stiff spring.

#3: For small molars, flat beaks.

#W3: Similar to #3, but wingless.

#4: Similar to #3, but with festooned beaks.

#5: For lower molars.

#56: Similar to #5, but with festooned beaks.

#56A: For large lowre molars, similar to #56, but with larger beaks.

#7: For lower molars with flat beaks.

#7A: For very big lower molars.

#26



#28



#30



#31



#51



#138



#139



#200



#201



#202



#203



#204



#205



#3



#W3



#4



#5



#56



#56A



#7



#7A



6 MOLAR CLAMPS

Autoclave 132°C

#7B: For very big molars

#8: For lower molars with round beaks.

#8A: For lower roots with strong festooned beaks.

#W8A: Similar to #8A, but wingless.

#10: For small lower molars, right and left.

#11: For small lower molars, right and left.

#12A: For third molars, serrated beaks.

#13A: For third molars, serrated beaks.

#14: For half erupted molars.

#14A: Similar to #14, but with larger beaks.

#14AN: Similar to #14A, but rather flat beaks.

#15: For big lower molars, wide beak strong spring.

#18: For upper and lower molars, round type beak.

#24: For deep buccal cavities of molars, wider fringe of buccal side.

7 ANTERIOR CLAMP

Autoclave 132°C

#210: For centrals.

#211: Universal for labial cavities on anterior.








#212: Cervical clamp for labial cavities.

#214: Cervical clamp

#6: For labial cavities.

#9: Most popular anterior clamp

#90N: Strong anterior clamps.



8 BICUSPID CLAMPS

Autoclave 132°C

#2 : For lower bicuspid with flat beaks

#2A : Similar to #2 with larger beaks.

#W2A : Similar to #2A without wings.

#0 : For small bicuspid
pedo central.

#00 : Smaller than #0

#22: Similar to #207 but wingless

#27: Similar to #206 but wingless
festooned beaks.

#29: For upper and lower bicuspid
broad beaks for wider sight.

#206: For upper and lower
bicuspid with festooned
beaks.

#207: For upper and lower
bicuspid with flat beaks.

#208: For big bicuspid with
similar pattern of #207

#209: For lower bicuspid with flat
beaks.

#1: For roots, deep festooned
beaks.

#1A: For cervical cavities of
bicuspid.

9 SPECIAL CLAMP

Autoclave 132°C

#P-1: For child first molars.

#P-2: For child first molars.

#U-67: For upper second molars.

#G-1: Not for rubberdam, used
with gauze or cotton
rolls.

#G-2: Very big wings is for tongue
and cheek press.

#2



#2A



#W2A



#0



#00



#22



#27



#29



#206



#207



#208



#209



#1



#1A



#P-1



#P-2



#U-67



#G-1



#G-2



10 CLAMP BOARD

Autoclave 132°C

*You can choose this board with or without clamps.

#8 : 90mm x 90mm

With 8 clamps :

#9 x 2, #2A, #W2A,
#5, #56, #203, #204

#12 : 150mm x 90mm

With 12 clamps :

#210, #211, #0, #2, #2A,
#206, #207, #5, #56, #7,
#201, #202

#18 : 160mm x 120mm

With 18 clamps :

#22, #27, #29, #208, #209,
#18, #138, #139, #200, #201
#202, #203, #204, #205.

I - SERIES : 260mm x 165mm

With 28 clamps.

#6, #9, #90N, #00, #0, #1
#1A, #2, #2A, #W2A, #3,
#W3, #4, #5, #7, #7A, #7B
#8, #8A, #W8A, #10, #11
#12A, #13A, #14, #14A,
#56, #56A

S - SERIES : 260mm x 165mm


With 31 clamps.

#210, #211, #212, #15, #18
#22, #24, #25, #26, #27, #28
#29, #30, #31, #51, #138
#139, #200, #201, #202
#203, #204, #205, #206
#207, #208, #209,
#G-1, #G-2, #P-1, #P-2

#8




#12



#18



I - SERIES



S - SERIES

